


# Custom OCR using YOLOv3 & Tesseract



# Project Overview

01 Combines YOLOv3 and  
Tesseract OCR

02 Detects and extracts structured  
text from lab reports

03 Generates downloadable CSV  
and Excel outputs

04 Streamlit-based user interface

05 Cloud deployment with AWS &  
Google Colab






- **YOLOv3**: Detects structured fields in lab reports
- **Tesseract OCR**: Extracts readable text from detected regions
- **OpenCV**: Helps with image preprocessing
- **Streamlit**: Offers a user-friendly web interface
- **Pandas & ExcelWriter**: Format and export clean data

## Tech Stack & Dependencies

### Core Tools Used



# Workflow Architecture





- Convert image to grayscale for contrast
- Apply Gaussian blur and thresholding
- Use bitwise inversion to standardize text
- Clean input regions before Tesseract OCR
- Improves clarity and output structure

## OCR & Image Preprocessing

**Enhancing  
Recognition  
Accuracy**





- Upload lab report images via browser
- Real-time preview of detection output
- View extracted text and patient data
- Download structured data as CSV or Excel
- Accessible on desktop or mobile

## Streamlit App Interface

**User-Friendly  
Experience**





# Lab Report Upload Interface

## Streamlit UI Upload

Users can upload lab report images directly through the Streamlit interface. The system supports JPG, JPEG, and PNG formats.

## Live Preview

The uploaded image is previewed on the screen before processing begins. This helps ensure the correct file is selected.

## File Management

Streamlit handles file input cleanly and gives the user clear visual feedback for the file upload and status.





# OCR Output Table

## Structured Data

After processing, extracted lab data is displayed in a structured table showing test names, values, units, and reference ranges.

## Patient Information

Patient name is extracted separately and highlighted above the results for clarity.

## Export Options

Users can download the OCR results as CSV or Excel with a single click, making data portable and editable.

Patient Name

ASHMIK SENGUPTA (SYM)

OCR Completed

TEST NAME	VALUE	UNIT	REFERENCE
1. TSH	0.5	mIU/L	0.1-0.5
2. FT4	1.8	ng/dL	0.8-1.2
3. T4	1.8	ng/dL	0.8-1.2
4. T4	1.8	ng/dL	0.8-1.2

Download CSV

Download Excel






YOLOv3 identifies and labels key fields in the report such as patient name, test name, values, and units with bounding boxes.

Each label is accompanied by a confidence score to evaluate detection accuracy.

These regions are then cropped and passed to Tesseract for final text recognition.





- ☐ Colab for development & GPU-based training
- ☐ Streamlit Cloud for hosting the UI
- ☐ Results downloadable in real-time
- 📦 Future: Dockerize and deploy on AWS EC2 or GCP VM
- 🔒 Can be scaled to private/secure environments

## Cloud Deployment Options

Bringing It to Production





- ☐ Support for multi-page PDF parsing
- 🛠️ Handwritten text recognition using deep learning models
- 🌐 Deploy as REST API for broader integration
- ✅ Compare accuracy with AWS Textract and Google Vision
- 🌍 Add multilingual OCR support for global use cases

## Future Work & Expansion

### Opportunities for Enhancement





Custom OCR Demo

# Thank You!



Any Questions?

Thank you for your time and attention.  
I'm happy to walk you through the  
demo, code, or discuss future  
improvements.

